


Anti-Otub2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Dc9

Item Number: CM0020495



Introduction

Anti-OTUB2 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human ubiquitin thioesterase OTUB2, a deubiquitinase involved in protein turnover. Ideal for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	OTUB2 Rabbit pAb
Remarks/Alias	OTB2, OTU2, C14orf137, OTUB2
Species	Human
GeneID	78990
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-73 of human OTUB2 (NP_075601.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96DC9
Protein Weight	27kDa
Shipping	Ice bag