

Anti-Otud4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q01804

Item Number: CM0023797

Introduction

Rabbit polyclonal antibody targeting OTUD4 (HIN1), validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Investigate OTUD4's role in deubiquitination, inflammatory signaling, and DNA repair.

[Learn More](#)

Parameter	Specification
Product Name	OTUD4 Rabbit pAb
Remarks/Alias	HIN1; DUBA6; HSHIN1; OTUD4
Species	Human
GeneID (Human)	54726
GeneID	54726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-500 of human OTUD4 (NP_001096123.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01804
Protein Weight	124kDa
Shipping	Ice bag