

Anti-Otud7B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6Gqq9

Item Number: CM0021198



Introduction

OTUD7B rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat OTUD7B. Ideal for NF-kB pathway and inflammation research.

[Learn More](#)

Parameter	Specification
Product Name	OTUD7B Rabbit pAb
Remarks/Alias	ZA20D1; CEZANNE; OTUD7B
Species	Human
GeneID (Human)	56957
GeneID	56957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-700 of human OTUD7B (NP_064590.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6GQQ9
Protein Weight	93kDa
Shipping	Ice bag