

Anti-Otx1 Monoclonal Antibody For Wb And Elisa - P32242

Item Number: CM0003298

Introduction

High-quality rabbit monoclonal antibody targeting human OTX1 (Homeobox protein OTX1). Suitable for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for brain development studies.

[Learn More](#)

Parameter	Value
Product Name	OTX1 Rabbit mAb
Species	Human
GeneID	5013
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human OTX1 (NP_055377.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P32242
Protein Weight	37kDa
Shipping	Ice bag