

Anti-Otx2 Rabbit Monoclonal Antibody For Wb, Elisa - P32243

Item Number: CM0011794

Introduction

High-quality OTX2 rabbit monoclonal antibody for western blot and ELISA applications, validated for human, mouse, and rat. Supports studies in brain development and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	OTX2 Rabbit mAb
Remarks/Alias	CPHD6; MCOP55; OTX2
Species	Human
GeneID (Human)	5015
GeneID	5015
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 190-289 of human OTX2 (P32243)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P32243
Protein Weight	32kDa
Shipping	Ice bag