

Anti-Otx2 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P32243

Item Number: CM0017144

Introduction

High-quality anti-OTX2 rabbit polyclonal antibody for immunohistochemistry, immunofluorescence, and ELISA. Cross-reacts with mouse and rat OTX2. Ideal for studying brain and sense organ development.

[Learn More](#)

Parameter	Value
Product Name	OTX2 Rabbit pAb
Remarks/Alias	CPHD6; MCOP55; OTX2
Species	Human
Gene ID (Human)	5015
Gene ID	5015
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 140-280 of human OTX2 (NP_068374.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P32243
Protein Weight	32kDa
Shipping	Ice bag