

Anti-Ovgp1 Rabbit Polyclonal Antibody For Wb And Elisa - Q12889

Item Number: CM0013486

Introduction

Rabbit polyclonal antibody targeting human OVGP1, validated for Western blot and ELISA in human, mouse, and rat. Ideal for research on fertilization, early embryonic development, and oviductal biology.

[Learn More](#)

Parameter	Value
Product Name	OVGP1 Rabbit pAb
Remarks/Alias	EGP, OGP, MUC9, CHIT5, OVGP1
Species	Human
GeneID (Human)	5016
GeneID	5016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-260 of human OVGP1 (NP_002548.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12889
Protein Weight	75kDa
Shipping	Ice bag