

Anti-Ovol2 Polyclonal Antibody For Wb, Elisa - Q9Brp0

Item Number: CM0020023

Introduction

Rabbit polyclonal antibody targeting human OVOL2 (ZNF339). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for epithelial-mesenchymal transition, developmental biology, and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	OVOL2 Rabbit pAb
Remarks/Alias	CHED; CHED1; CHED2; PPCD1; ZNF339; EUROIMAGE566589; OVOL2
Species	Human
Gene ID (Human)	58495
Gene ID	58495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-98 of human OVOL2 (NP_067043.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BRP0
Protein Weight	30kDa
Shipping	Ice bag