


Anti-Ox40L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P23510

Item Number: CM0021687



Introduction

Rabbit polyclonal antibody against human OX40L (CD252, TNFSF4), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Suitable for T-cell co-stimulation and cytokine studies.

[Learn More](#)

Parameter	Specification
Product Name	OX40L Rabbit pAb
Remarks/Alias	GP34; CD252; OX40L; TXGP1; CD134L; OX-40L; TNLG2B; OX40L
Species	Human
Gene ID (Human)	7292
Gene ID	7292
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 51-183 of human OX40L (NP_003317.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23510
Protein Weight	21kDa
Shipping	Ice bag