

Anti-Oxct1 Rabbit Polyclonal Antibody For Wb, Elisa - P55809

Item Number: CM0030265



Introduction

Rabbit polyclonal antibody against OXCT1 (SCOT), a key enzyme for ketone body catabolism, validated for WB and ELISA in human, mouse, and rat samples. Suitable for metabolic research and mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	OXCT1 Rabbit pAb
Remarks/Alias	OXCT; SCOT; OXCT1
Species	Human
Gene ID	5019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 261-520 of human OXCT1 (NP_000427.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55809
Protein Weight	56kDa
Shipping	Ice bag