

Anti-Oxct2 Polyclonal Antibody For Wb, If-P, Elisa - Q9Byc2

Item Number: CM0018810

Introduction

Rabbit polyclonal antibody against human OXCT2, a mitochondrial ketone body catabolism enzyme. Validated for WB, IF-P, and ELISA applications with cross-reactivity to human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	OXCT2 Rabbit pAb
Remarks/Alias	SCOTT; FKSG25; OXCT2
Species	Human
GeneID (Human)	64064
GeneID	64064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-510 of human OXCT2 (NP_071403.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYC2
Protein Weight	56kDa
Shipping	Ice bag