

Anti-Oxr1 Monoclonal Antibody For Wb, If-P, Elisa - Q8N573

Item Number: CM0006535

Introduction

Recombinant anti-OXR1 rabbit monoclonal antibody validated for WB, IF-P, ELISA. Detects 24-98 kDa isoforms in human, mouse, rat. Ideal for studying oxidative stress, V-ATPase regulation, and neuroprotection.

[Learn More](#)

Parameter	Value
Product Name	OXR1 Rabbit mAb
Remarks/Alias	TLDC3; CHEGDD; Nbla00307; OXR1
Species	Human
GeneID (Human)	55074
GeneID	55074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 705-874 of human OXR1 (NP_001185461.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N573
Protein Weight	24-98 kDa
Shipping	Ice bag