

Anti-Oxr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N573

Item Number: CM0019081

Introduction

Rabbit polyclonal antibody against human OXR1 (Oxidation resistance protein 1), validated for Western blot and ELISA in human, mouse, and rat samples. Useful for studying oxidative stress, Golgi pH regulation, and neuroprotection.

[Learn More](#)

Parameter	Value
Product Name	OXR1 Rabbit pAb
Remarks/Alias	TLDC3; CHEGDD; Nbla00307; OXR1
Species	Human
GeneID (Human)	55074
GeneID	55074
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 705-874 of human OXR1. (NP_001185461.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N573
Protein Weight	98kDa
Shipping	Ice bag