


Anti-Oxsm Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nwu1

Item Number: CM0020621



Introduction

Validate OXSM rabbit polyclonal antibody for Western blot and ELISA. Reacts with human, mouse, rat 3-oxoacyl-ACP synthase, mitochondrial. Supports fatty acid biosynthesis, lipid metabolism, and mitochondrial function studies.

[Learn More](#)

Parameter	Specification
Product Name	OXSM Rabbit pAb
Remarks/Alias	KS; CEM1; KASI; FASN2D; OXSM
Species	Human
Gene ID (Human)	54995
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-459 of human OXSM (NP_060367.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWU1
Protein Weight	49kDa
Shipping	Ice bag