


Anti-Oxsr1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95747

Item Number: CM0018351



Introduction

Rabbit polyclonal antibody targeting OXSR1 (OSR1), a serine/threonine kinase in the WNK-OSR1 pathway. Verified for WB, IF/ICC, ELISA in human, mouse, and rat. Suitable for ion transport, stress response, and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	OXSR1 Rabbit pAb
Remarks/Alias	OSR1; OXSR1
Species	Human
Gene ID (Human)	9943
Gene ID	9943
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OXSR1 (NP_005100.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95747
Protein Weight	58kDa
Shipping	Ice bag