

Anti-Oxt Rabbit Polyclonal Antibody For Wb, Elisa - P01178

Item Number: CM0016414



Introduction

High-quality rabbit polyclonal antibody against human OXT (oxytocin-neurophysin 1 proprotein). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for social behavior and reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	OXT Rabbit pAb
Remarks/Alias	OT; OT-NPI; OXT-NPI; OXT
Species	Human
Gene ID	5020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human OXT (NP_000906.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01178
Protein Weight	13kDa
Shipping	Ice bag