


Anti-Oxytocin Receptor Monoclonal Antibody For Wb, Elisa - P30559

Item Number: CM0003653



Introduction

Rabbit monoclonal antibody against human oxytocin receptor (OXTR, P30559). For WB and ELISA applications, with cross-reactivity to mouse and rat. Ideal for research on neuropeptide hormone signaling, social behavior, and receptor pharmacology.

[Learn More](#)

Parameter	Value
Product Name	Oxytocin Receptor Rabbit mAb
Remarks/Alias	OTR; OT-R; Oxytocin Receptor
Species	Human
Gene ID (Human)	5021
Gene ID	5021
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human Oxytocin Receptor (P30559)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30559
Protein Weight	43kDa
Shipping	Ice bag