


## Anti-P-Cadherin Rabbit Monoclonal Antibody For Wb, If, Fc And Elisa - P22223

Item Number: CM0004106



### Introduction

High-quality rabbit monoclonal antibody against human P-cadherin (CDH3), validated for WB, IF/ICC, IF-P, FC, ELISA. Recognizes endogenous P-cadherin protein (~87/91 kDa). Ideal for cell adhesion and carcinoma research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | P-cadherin Rabbit mAb                                                                                                   |
| Remarks/Alias    | CDHP; HJMD; PCAD                                                                                                        |
| Species          | Human                                                                                                                   |
| Gene ID (Human)  | 1001                                                                                                                    |
| Gene ID          | 1001                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 108-654 of human P-cadherin (NP_001784.2) |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, IF-P, FC, ELISA                                                                                             |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | P22223                                                                                                                  |
| Protein Weight   | 87kDa/91kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                                 |