


Anti-P Glycoprotein Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P08183

Item Number: CM0005227



Introduction

Rabbit monoclonal antibody targeting P Glycoprotein (ABCB1/MDR1) with validated applications in WB, IF-P, IHC-P, and ELISA. Cross-reactive with human, mouse, and rat. Suitable for multidrug resistance and transporter studies.

[Learn More](#)

Parameter	Value
Product Name	P Glycoprotein Rabbit mAb
Remarks/Alias	CLCS; MDR1; P-GP; PGY1; ABC20; CD243; GP170; p-170
Species	Human
GeneID (human)	5243
GeneID	5243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 603-692 of human P Glycoprotein (NP_000918.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08183
Protein Weight	141kDa
Shipping	Ice bag