

Anti-P Glycoprotein Polyclonal Antibody For Wb, If/Icc, Elisa - P08183

Item Number: CM0018597



Introduction

High-quality anti-P glycoprotein rabbit polyclonal antibody for WB, IF/ICC, and ELISA. Detects human, mouse, and rat ABCB1 (MDR1) protein. Suitable for multidrug resistance and drug transport research.

[Learn More](#)

Parameter	Value
Product Name	P Glycoprotein Rabbit pAb
Remarks/Alias	CLCS; MDR1; P-GP; PGY1; ABC20; CD243; GP170; p-170; P Glycoprotein
Species	Human
GeneID	5243
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 603-692 of human P Glycoprotein (NP_000918.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08183
Protein Weight	141kDa
Shipping	Ice bag