

Anti-P Glycoprotein Rabbit Monoclonal Antibody For Wb, Elisa - P08183

Item Number: CM0002814



Introduction

Rabbit monoclonal antibody against human P Glycoprotein (ABCB1/MDR1), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for multidrug resistance and drug transport studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] P Glycoprotein Rabbit mAb
Remarks/Alias	CLCS; MDR1; P-GP; PGY1; ABC20; CD243; GP170; p-170; P Glycoprotein
Species	Human
GeneID (Human)	5243
GeneID	5243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 393-710 of human P Glycoprotein (P08183).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08183
Protein Weight	141kDa
Shipping	Ice bag