

Anti-P Glycoprotein Rabbit Pab - P08183

Item Number: CM0018213



Introduction

Rabbit polyclonal antibody targeting human P Glycoprotein (ABCB1). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse. Key tool for multidrug resistance and ABC transporter studies.

[Learn More](#)

Parameter	Value
Product Name	P Glycoprotein Rabbit pAb
Remarks/Alias	CLCS; MDR1; P-GP; PGY1; ABC20; CD243; GP170; p-170; P Glycoprotein
Species	Human
GeneID (Human)	5243
GeneID	5243
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 590-689 of human P Glycoprotein (NP_000918.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08183
Protein Weight	141kDa
Shipping	Ice bag