

Anti-P150 Caf1 (Caf-1 P150 Subunit) Rabbit Polyclonal Antibody For Wb And Elisa - Q13111

Item Number: CM0022717



Introduction

Rabbit polyclonal antibody targeting human p150 CAF1 (CHAF1A), the p150 subunit of chromatin assembly factor-1. Validated for Western blot and ELISA. Suitable for chromatin assembly, DNA replication, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	p150 CAF1 Rabbit pAb
Remarks/Alias	CAF1; P150; CAF-1; CAF1B; CAF1P150; p150 CAF1
Species	Human
GeneID (human)	10036
GeneID	10036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human p150 CAF1 (NP_005474.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13111
Protein Weight	107kDa
Shipping	Ice bag