

Anti-P150 Caf1 Rabbit Monoclonal Antibody For Wb, Elisa - Q13111

Item Number: CM0009898



Introduction

High-quality rabbit monoclonal antibody against human p150 CAF1 (CHAF1A), validated for Western blot and ELISA. Ideal for studying chromatin assembly and DNA replication/repair processes.

[Learn More](#)

Parameter	Value
Product Name	p150 CAF1 Rabbit mAb
Remarks/Alias	CAF1; P150; CAF-1; CAF1B; CAF1P150; p150 CAF1
Species	Human
GeneID (Human)	10036
GeneID	10036
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13111
Protein Weight	107kDa
Shipping	Ice bag