


Anti-P16-Arc/Arpc5 Polyclonal Antibody For Wb, Elisa - O15511

Item Number: CM0029316



Introduction

Rabbit polyclonal antibody against human ARPC5 (p16-ARC), validated for WB and ELISA. Cross-reacts with mouse and rat. Essential tool for studying Arp2/3 complex, actin dynamics, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	p16-ARC/ARC16/ARPC5 Rabbit pAb
Remarks/Alias	ARC16; p16-Arc; dj127C7.3; p16-ARC/ARC16/ARPC5
Species	Human
GeneID	10092
Immunogen	Recombinant fusion protein containing amino acids 1-154 of human ARPC5 (NP_001257368.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15511
Protein Weight	16kDa
Shipping	Ice bag