

## Anti-P16 Ink4A Polyclonal Antibody For Wb - P51480

Item Number: CM0019164

### Introduction

High-quality rabbit polyclonal antibody targeting the N-terminus of p16 INK4A. Validated for WB and ELISA with cross-reactivity to human, mouse, rat. Ideal for studying cell cycle regulation and tumor suppression.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | p16 INK4A Rabbit pAb                                                                                                      |
| Remarks/Alias    | Arf; p16; MTS1; Pctr1; p19ARF; p16INK4a; p19<ARF>; ARF-INK4a; INK4a-ARF; Ink4a/Arf; p16(INK4a); p16 INK4A                 |
| Species          | Mouse                                                                                                                     |
| Gene ID (Human)  | 1029                                                                                                                      |
| Gene ID          | 12578                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of mouse p16 INK4A. (NP_001035744.1). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, ELISA                                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | P51480                                                                                                                    |
| Protein Weight   | 17KDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |