


## Anti-P19Arf Rabbit Polyclonal Antibody For Wb, Elisa - Q64364

Item Number: CM0021935



### Introduction

Rabbit polyclonal antibody against mouse p19ARF (Tumor suppressor ARF), validated for WB and ELISA. Cross-reacts with human. Suitable for cell cycle and tumor suppression studies.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | p19ARF Rabbit pAb                                                                                       |
| Remarks/Alias    | Arf; MTS1; Pctr1; p19ARF; p19<ARF>; ARF-INK4a; INK4a-ARF; Ink4a/Arf;                                    |
| Species          | Mouse                                                                                                   |
| GeneID (Human)   | 1029                                                                                                    |
| GeneID           | 12578                                                                                                   |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse p19ARF (NP_034007.1). |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, ELISA                                                                                               |
| Cross Reactivity | Human                                                                                                   |
| SWISS            | Q64364                                                                                                  |
| Protein Weight   | 19 kDa                                                                                                  |
| Shipping         | Ice bag                                                                                                 |