


Anti-P23/Ptges3 Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q15185

Item Number: CM0017453



Introduction

High-quality rabbit polyclonal antibody targeting P23/PTGES3 (cPGES), validated for IHC-P and ELISA. Detects human, mouse, and rat protein; key reagent for studying prostaglandin E2 synthesis and chaperone functions.

[Learn More](#)

Parameter	Value
Product Name	P23/PTGES3 Rabbit pAb
Remarks/Alias	P23, TEBP, cPGES
Species	Human
Gene ID (Human)	10728
Gene ID	10728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human P23/PTGES3 (NP_006592.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15185
Protein Weight	19kDa
Shipping	Ice bag