


Anti-P23/Ptges3 Rabbit Polyclonal Antibody For Wb And Elisa - Q15185

Item Number: CM0030506



Introduction

Rabbit polyclonal antibody targeting human PTGES3 (Prostaglandin E synthase 3, cPGES). Validated for Western blotting and ELISA; cross-reacts with mouse and rat. Ideal for studying prostaglandin biosynthesis, chaperone function, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	P23/PTGES3 Rabbit pAb
Remarks/Alias	P23; TEBP; cPGES; P23/PTGES3
Species	Human
Gene ID (Human)	10728
Gene ID	10728
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human P23/P23/PTGES3 (NP_006592.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15185
Protein Weight	19kDa
Shipping	Ice bag