

Anti-P2Rx3 Polyclonal Antibody For Wb, Elisa - P56373

Item Number: CM0017305



Introduction

Rabbit polyclonal antibody specific for human P2RX3 (P2X purinoceptor 3), suitable for WB and ELISA. Detects mouse and rat P2RX3. Target is an ATP-gated ion channel involved in sensory processing and pain perception.

[Learn More](#)

Parameter	Value
Product Name	P2RX3 Rabbit pAb
Remarks/Alias	P2X3; P2RX3
Species	Human
Gene ID (Human)	5024
Gene ID	5024
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P56373
Protein Weight	44kDa
Shipping	Ice bag