

Anti-P2Rx4 Polyclonal Antibody For If/Icc, If-P, Elisa - Q99571

Item Number: CM0024352



Introduction

Rabbit polyclonal antibody against human P2RX4 (P2X4 purinoceptor, ATP-gated cation channel). Validated in IF/ICC, IF-P, ELISA. Detects human and rat P2RX4. Ideal for studying immune response and neuroinflammation.

[Learn More](#)

Parameter	Value
Product Name	P2RX4 Rabbit pAb
Remarks/Alias	P2X4; P2X4R; P2RX4
Species	Human
GeneID	5025
Immunogen	Synthetic peptide: a synthetic peptide corresponding to a sequence within amino acids 111-211 of human P2RX4 (NP_002551.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q99571
Protein Weight	43kDa
Shipping	Ice bag