

Anti-P2Rx4 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q99571

Item Number: CM0004889



Introduction

High-quality rabbit monoclonal antibody against human P2RX4 (P2X4 purinoceptor). Validated for WB, IF-P, ELISA; cross-reacts with human, mouse, rat. A useful tool for studying ATP-gated ion channels in immunity and neuroscience.

[Learn More](#)

Parameter	Value
Product Name	P2RX4 Rabbit mAb
Remarks/Alias	P2X4; P2X4R; P2RX4
Species	Human
GeneID	5025
Immunogen	Recombinant fusion protein containing amino acids 55-338 of human P2RX4 (NP_002551.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99571
Protein Weight	43kDa
Shipping	Ice bag