

Anti-P2Rx4 Rabbit Polymab® For Wb, If/Icc, If-P, Ihc-P, Elisa - Q99571

Item Number: CM0011340



Introduction

High-sensitivity rabbit recombinant monoclonal anti-P2RX4 antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Detects human, mouse, rat P2X4 purinoceptor; ideal for immunology and neuroscience research. ~43 kDa.

[Learn More](#)

Parameter	Value
Product Name	P2RX4 Rabbit PolymAb®
Remarks/Alias	P2X4; P2X4R; P2RX4
Species	Human
GeneID (Human)	5025
GeneID	5025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-338 of human P2RX4 (NP_110378.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99571
Protein Weight	43kDa
Shipping	Ice bag