

Anti-P2Rx5 Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q93086

Item Number: CM0016007

Introduction

Rabbit polyclonal antibody targeting human P2RX5 (P2X5, LRH-1), validated for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying ATP-gated ion channels in brain and immune tissues.

[Learn More](#)

Parameter	Value
Product Name	P2RX5 Rabbit pAb
Remarks/Alias	P2X5; LRH-1; P2X5R; P2RX5
Species	Human
GeneID (Human)	5026
GeneID	5026
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-422 of human P2RX5 (NP_002552.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93086
Protein Weight	47kDa
Shipping	Ice bag