

Anti-P2Rx7 Polyclonal Antibody For Wb, Elisa - Q99572

Item Number: CM0015480



Introduction

P2RX7 Rabbit pAb is a polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Target: P2X7 receptor (Q99572), an ATP-gated ion channel involved in inflammation and cell death.

[Learn More](#)

Parameter	Value
Product Name	P2RX7 Rabbit pAb
Remarks/Alias	P2X7; P2RX7
Species	Human
Gene ID (Human)	5027
Gene ID	5027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-334 of human P2RX7 (NP_002553.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99572
Protein Weight	69kDa
Shipping	Ice bag