

Anti-P2Ry1 Rabbit Polyclonal Antibody For Wb And Elisa - P47900

Item Number: CM0029089



Introduction

Rabbit polyclonal antibody against human P2RY1 (P2Y1 purinoceptor 1), validated for WB and ELISA. Reacts with human, mouse, and rat. Useful for studying ADP-mediated platelet aggregation and hemostasis.

[Learn More](#)

Product Name	Remarks/Alias	Species	GeneID (Human)	GeneID	Immunogen	Source	Category	Application	Cross Reactivity	SWISS	Protein Weight	Shipping
P2RY1 Rabbit pAb	P2Y1, SARCC, P2RY1	Human	5028	5028	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 274-373 of human P2RY1 (NP_002554.1).	Rabbit	Polyclonal Antibodies	WB, ELISA	Human, Mouse, Rat	P47900	42kDa	Ice bag