

Anti-P2Ry10 Rabbit Polyclonal Antibody For Wb, Elisa - 000398

Item Number: CM0024036



Introduction

Rabbit polyclonal antibody against P2RY10 (putative P2Y purinoceptor 10), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Useful for studying G protein-coupled purinergic signaling.

[Learn More](#)

Parameter	Specification
Product Name	P2RY10 Rabbit pAb
Remarks/Alias	LPS2; P2Y10; LYPSR2; P2RY10
Species	Human
GeneID (Human)	27334
GeneID	27334
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human P2RY10 (NP_001311147.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00398
Protein Weight	39kDa
Shipping	Ice bag