

Anti-P2Ry11 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q96G91

Item Number: CM0023120

Introduction

High-quality anti-P2RY11 rabbit polyclonal antibody for human and mouse samples. Validated for IF/ICC and ELISA applications. Detects the P2Y purinoceptor 11 (P2Y11) receptor, a G protein-coupled receptor, at approximately 40 kDa.

[Learn More](#)

Parameter	Value
Product Name	P2RY11 Rabbit pAb
Remarks/Alias	P2RY11; P2Y11; P2Y purinoceptor 11
Species	Human
Gene ID (human)	5032
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human P2RY11 (NP_002557.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96G91
Protein Weight	40kDa
Shipping	Ice bag