

Anti-P2Ry11 Rabbit Polyclonal Antibody For Wb, Elisa - Q96G91

Item Number: CM0023711

Introduction

Rabbit polyclonal antibody targeting human P2Y purinoceptor 11 (P2RY11), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for G protein-coupled receptor signaling research.

[Learn More](#)

Parameter	Value
Product Name	P2RY11 Rabbit pAb
Remarks/Alias	P2Y11; P2RY11
Species	Human
GeneID (human)	5032
GeneID	5032
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 300-379 of human P2RY11 (NP_002557.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96G91
Protein Weight	40kDa
Shipping	Ice bag