


Anti-P2Ry12 Polyclonal Antibody For Wb, If-P, Elisa - Q9H244

Item Number: CM0022821



Introduction

Rabbit polyclonal antibody targeting human P2RY12 (P2Y purinoceptor 12), validated for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Essential for platelet aggregation and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	P2RY12 Rabbit pAb
Remarks/Alias	HORK3; P2Y12; ADPG-R; BDPLT8; SP1999; P2T(AC); P2Y(AC); P2Y(12)R; P2Y(ADP); P2Y(cyc); P2RY12
Species	Human
GeneID	64805
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 243-342 of human P2RY12 (NP_795345.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H244
Protein Weight	39kDa
Shipping	Ice bag