

Anti-P2Ry14 Rabbit Polyclonal Antibody For Wb, Ihc-P, And Elisa - Q15391

Item Number: CM0027869



Introduction

Rabbit polyclonal antibody against human P2RY14 (P2Y14 receptor, GPR105), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying UDP-glucose signaling and GPCR biology.

[Learn More](#)

Parameter	Value
Product Name	P2RY14 Rabbit pAb
Remarks/Alias	P2Y14; BPR105; GPR105; P2RY14
Species	Human
GeneID (human)	9934
GeneID	9934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-338 of human P2RY14 (NP_055694.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15391
Protein Weight	39kDa
Shipping	Ice bag