


Anti-P2Ry2 Polyclonal Antibody For Wb, Elisa - P41231

Item Number: CM0013723



Introduction

Rabbit polyclonal antibody targeting human P2RY2 (P41231), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	P2RY2 Rabbit pAb
Remarks/Alias	P2U; HP2U; P2U1; P2UR; P2Y2; P2RU1; P2Y2R; P2RY2
Species	Human
GeneID (human)	5029
GeneID	5029
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 278-377 of human P2RY2 (NP_788086.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41231
Protein Weight	42kDa
Shipping	Ice bag