


Anti-P2Ry2 Rabbit Polyclonal Antibody For Wb, Elisa - P41231

Item Number: CM0030382



Introduction

Rabbit polyclonal antibody targeting human P2Y purinoceptor 2 (P2RY2), validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for GPCR signaling studies. Approx. 42 kDa.

[Learn More](#)

Parameter	Value
Product Name	P2RY2 Rabbit pAb
Remarks/Alias	P2U; HP2U; P2U1; P2UR; P2Y2; P2RU1; P2Y2R; P2RY2
Species	Human
GeneID (human)	5029
GeneID	5029
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P41231
Protein Weight	42kDa
Shipping	Ice bag