


Anti-P38 Mapk (Ko Validated) Rabbit Monoclonal Antibody For Wb, Elisa - Q16539

Item Number: CM0008646



Introduction

KO validated anti-p38 MAPK rabbit monoclonal antibody for Western blot and ELISA. Detects human, mouse, rat p38 MAPK (41 kDa). Suitable for MAPK signaling, inflammation, and stress response research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] p38 MAPK Rabbit mAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; PK
Species	Human
GeneID (Human)	1432
GeneID	1432
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag