


Anti-P38 Mapk Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q16539

Item Number: CM0013074



Introduction

Rabbit polyclonal antibody against human p38 MAPK (MAPK14), suitable for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, and rat. Ideal for stress signaling, inflammation, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	p38 MAPK Rabbit pAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA
Species	Human
Gene ID	1432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 261-360 of human p38 MAPK (NP_620581.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag