


Anti-P38 Mapk Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q16539

Item Number: CM0005186



Introduction

Anti-p38 MAPK Rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat MAPK14 protein. Suitable for studying stress-activated MAP kinase cascades.

[Learn More](#)

Parameter	Value
Product Name	p38 MAPK Rabbit mAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; p38 MAPK
Species	Human
Gene ID	1432
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 100-200 of human p38 MAPK (Q16539)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag