


Anti-P38 Mapk Rabbit Polyclonal Antibody For Wb - Q16539

Item Number: CM0016520



Introduction

High-quality rabbit polyclonal antibody against human p38 MAPK (SWISS: Q16539). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Detects a 41 kDa protein. Suitable for inflammation, apoptosis, and cellular stress research.

[Learn More](#)

Parameter	Specification
Product Name	p38 MAPK Rabbit pAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; p38 MAPK
Species	Human
GeneID (Human)	1432
GeneID	1432
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 261-360 of human p38 MAPK (NP_620581.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag