

Anti-P47Phox/Ncf1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P14598

Item Number: CM0016184



Introduction

Rabbit polyclonal antibody against human p47phox/NCF1, validated for Western blot, immunofluorescence/immunocytochemistry, immunoprecipitation, and ELISA. Cross-reacts with mouse and rat. Useful for studying NADPH oxidase function and chronic granulomatous disease.

[Learn More](#)

Parameter	Value
Product Name	p47phox/NCF1 Rabbit pAb
Remarks/Alias	CGD1; NCF1A; NOXO2; p47phox; SH3PXD1A; p47phox/NCF1
Species	Human
GeneID	653361
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-390 of human p47phox/NCF1 (NP_000256.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14598
Protein Weight	45kDa
Shipping	Ice bag