

Anti-P47Phox/Ncf1 Rabbit Monoclonal Antibody For Wb, Elisa - P14598

Item Number: CM0007735



Introduction

Rabbit monoclonal antibody to human p47phox/NCF1 (P14598). Validated for WB and ELISA, detecting a 45 kDa protein essential for NADPH oxidase activity. Suitable for CGD studies.

[Learn More](#)

Parameter	Value
Product Name	p47phox/NCF1 Rabbit mAb
Remarks/Alias	CGD1; NCF1A; NOXO2; p47phox; SH3PXD1A; p47phox/NCF1
Species	Human
GeneID (Human)	653361
GeneID	653361
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 291-390 of human p47phox/NCF1 (P14598).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14598
Protein Weight	45kDa
Shipping	Ice bag