

Anti-P4Ha2 Rabbit Monoclonal Antibody For Wb, Elisa - O15460

Item Number: CM0006882



Introduction

Reliable rabbit monoclonal antibody targeting human P4HA2, validated for Western blot and ELISA. Suitable for research on collagen hydroxylation, fibrosis, and prolyl 4-hydroxylase activity.

[Learn More](#)

Parameter	Value
Product Name	P4HA2 Rabbit mAb
Remarks/Alias	MYP25; lncRNA-PE; P4HA2
Species	Human
GeneID (Human)	8974
GeneID	8974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 157-333 of human P4HA2 (NP_004190.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15460
Protein Weight	61 kDa
Shipping	Ice bag